English	Manuel d'utilisation	Version: A1
Disposable EEG cable		
User Manual	Version: A1	

1. Indication for Use/Intended Use

The Disposable EEG cable is intended for non-invasive use with recording and monitoring equipment, (active and reference), of Electroencephalograph (EEG), electromyography (EMG), and Evoked Potentials (EP).

2. CAUTION

Law restricts this device to sale by or on the order of a physician.

3. Device description

Item	P01-001-00566	P01-001-00567	P01-001-00568	P01-001-00569
Classification	Cable	Cable	Cable	Cable
Length (in/mm)	48/1220	72/1830	96/2440	120/3050
Lead number	10	14	14	14
Distal connector design	Cup	Cup	Cup	Cup
Minimum Impedance (Ω)	0.5415	0.8075	1.0925	1.349
Specified Impedance (Ω)	0.57	0.85	1.15	1.42
Maximum Impedance (Ω)	0.5985	0.8925	1.2075	1.491

Performance specifications:

Electrodes contained within a set have been verified to have low, balanced impedance (see table). As a result, the electrodes will contribute no more than 0.15% distortion of the measured voltage.

4. Operation Environment

- Environment temperature: 0°C ~40°C
- Relative humidity: 15% ~85%
- Atmospheric pressure: 86kpa ~106kpa

5. Storage

Storage conditions are as follows:

- Temperature: -10 to +40°C
- Relative humidity: 15 to 85%
- Atmospheric pressure: 86kpa ~106kpa

6. Safety & Performance Compliance

Safety: IEC 60601-1: 2012

Biocompatibility: ISO 10993-5: 2009, ISO 10993-10: 2010

7. Installation for use:

- Clean the patient’s skin where the signal is to be collected.
- After applying conductive paste, the product is fixed at the place where the signal is to be collected. Use only conductive paste that has been cleared by the relevant authorities in the user’s country.
- Connect the plug to the Electroencephalograph (EEG)/electromyography (EMG)/ Evoked Potentials (EP) device.
- Proper installation can be verified by reviewing recorded signals.

Warning: MR-unsafe!

- Do not expose the device to a magnetic resonance (MR) environment.
- Thermal injury and burns may occur due to the metal components of the device that can heat during MR scanning.
- The device may not function properly due to the strong magnetic and radiofrequency fields generated by the MR scanner.

8. Replacement and disposal

- Replacement: Check the appearance before use. If the cable/ jacket damaged when opening the packing bag, please replace with a new EEG cable.
- Remove the cable after the test is completed, assist the patient to clean the skin.
- The device should be disposed after use. Do not reuse the disposable device.
- Disposal of the cables should comply with local regulation.

9. Warning

- This cable is for use only with legal marketed host device.
- Do not use the device on the broken or damaged skin.
- Do not use the device for the purpose of electrical stimulation
- Do not use the device over 24 hours.
- Please always remember to read the user manual for the compatible device before using the EEG cable.
- Do not allow any liquid to get inside the various connectors on the interface cable.
- Do not reuse the disposable EEG cables.
- Do not use a cable or part of a cable if there is any risk to the patient (e.g. damaged insulating material).
- This disposable product is for single use only and is not intended to be disinfected or cleaned after patient use. Dispose of the product after patient use. Multi-patient use without disinfection between patients increases the risk of infection. This disposable product is not designed for reuse. Changke cannot guarantee the performance specifications once this product has been disinfected or cleaned after patient use.
- If, in relation to the use of the cable, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country.

10. Symbol explanation

	Caution		Manufacturer		Date of manufacture
	Type BF applied part	REF	Reference number	LOT	Batch code
	Do not re-use		Consult instructions for use		The device should be sent to the special agencies according to local regulations for separate collection after their useful life.
EC REP	EU representative	MR	MR-unsafe	Rx only (U.S.)	Federal (U.S.) Law restricts this device to sale by or on the order of a physician
CE	CE marking	MD	Medical device indicator		Biodegradable materials
#	Part number	P/N:	Part number		Use-by date
UDI	Unique Device Identification				

 Shenzhen Changke Connect Electronics Co., Ltd. Address: A1-2nd and A2-4th floor of factory building, No.87 of Hengping Road, Baoan Community, Yuanshan Street, Longgang District, Shenzhen, China

EC REP **CE** Wellkang Ltd. Enterprise Hub, NW Business Complex,1 Beraghmore Rd., Derry, BT48 8SE, Northern Ireland, UK Exclusively marketed by Philips

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